iM3® iM3 SP6 Ultrasonic Scaler & Micromotor

INSTRUCTIONS FOR USE AND
TECHNICAL MANUAL

ENGLISH

FRANCAIS

DEUTSCH

ESPAÑOL

iM3 SP6 Ultrasonic Scaler &
Micromotor

FOR VETERINARY USE ONLY

ART-SP6 Scaler /Polisher Combo Unit

Instruction Manual

Bonart Co., Ltd.

4F-11, No.3 Wuquan 1st Rd., Xinzhuang
Dist., New Taipei City 24892, Taiwan, R.O.C.
Tel: 886-2-22983980
Fax:886-2-22983981
Email: sales@bonartmed.com

COPYRIGHT© 2014 BONART CO., LTD.

BONART Co., Ltd. (BONART) owns all rights to this publication. This publication shall be used solely as a reference for operation, maintenance, and repair on BONART equipment. This document may not be reproduced or disseminated in any form, by any means, for any other purpose.

In the event of inadvertent or deliberate publication, BONART intends to enforce its rights to this manual under the copyright law. Copying materials presented in this manual is strictly prohibited without written permission by BONART Co., Ltd.

Product specifications and information presented in this manual are subject to change without prior notice.

**PROPERTY OF BONART CO., LTD.
ALL RIGHTS RESERVED**

ART is a registered trademark of BONART Co., Ltd. & its affiliates.

Document No. RD-02-14-06 Rev.1.0 (04/14)

Printed in Taiwan, R.O.C.

Operator Safety

Please read the entire manual prior to operating the unit.

The ART-SP6 Scaler/Polisher Combo Unit should be operated, maintained and repaired by qualified and properly trained personnel.

- ***NOTE, CAUTION, AND WARNING STATEMENTS***

NOTE : Provides tips and advice

CAUTION : Provides correct operating or maintenance procedures

WARNING : Alerts users of danger that may cause severe injury when proper procedures are not followed

- ***SYMBOLS***

Complies with MDD 93/42/EEC

BF equipment

Grounding terminal

Attention! Please read instructions

Warnings

Important!

- **Only trained and qualified doctors, dentists, and hygienists should operate the ART-SP6**
- **Patients with pacemakers cannot be treated with the ART-SP6**
- **Do not immerse the ART-SP6 in water!**
- **If the patient or operator is pregnant, consult a physician prior to performing dental scaling for safety purposes.**
- **The ART-SP6 requires a shielded AC power cord.**
- **Plug the power cord into a well-grounded outlet.**

Attention to users with cardiac pacemakers!


Patients with pacemakers should avoid treatment with the ART-SP6. Electronic appliances including razors, hair dryers, microwave ovens, TV receptors, and other medical equipment, such as the ART-SP6, may interfere with the performance of pacemakers.

For more information on this subject, please refer to the following articles:

- "Advances in Cardiac Pacemaker", The New York Academy of Sciences, Vol. 167, Article 2, pp. 515-1075
- "Electromagnetic Radiation Interference with Cardiac Pacemaker", U. S. Department of Health, Education and Welfare
- "The Individual with a Pacemaker in the Dental Environment", Journal of the American Dental Association, Vol. 91, No. 6, pp. 1224-1229

Please Read Before Operating:


- Make sure the device is connected to a grounded outlet. Failure to meet this requirement may severely damage the user and device.


- Place the device on a flat and stable surface. Placing the device on a tilted or unstable surface may degrade performance and also cause damage to the device.


- Do not modify the device. Modifications to the device will invalidate the safety codes and warranty while endangering the patient and operator.

- Do not place heavy objects on top of the power cord and keep the power cord away from intense heat.

- Unplug the device and call your local dealer or BONART for service and instructions if you observe any abnormalities while operating the device.

Preface

■ **About the ART-SP6 Scaler/Polisher Combo Unit**

The ART-SP6 is CE 0434 certified in compliance with the applicable requirements of the Council Directive 93/42/EEC. The reviewing council also classifies the ART-SP6 in accordance with EN60601-1 and EN60601-1-2 safety requirements.

The ART-SP6 utilizes piezoelectric ultrasonic waves generated by 4 ceramic piezoelectric plates located in the hand piece. These plates, subject to high frequency current begin to vibrate. The ART-SP6 contains a highly efficient transducer that provides an operating frequency of 29 KHz with minimal noise at a lower temperature, thus reducing the amount of cooling water needed during operation. The vibrations from the hand piece in combination with lower cooling water requirements produce the perfect device for removing tough calculus and tartar with ease.

■ **TECHNICAL SUPPORT**

For technical support and mail-in repairs in the U.S., please contact our U.S. division, Magpie Tech. Corp. at 1-(888) 526-6278, Monday through Friday, 9:00AM to 5:00PM (Pacific). For non US territories, please contact your local BONART representative.

■ **REPLACEMENT PARTS**

To order replacement parts, please contact an authorized dealer listed on our website, www.bonartmed.com

Table of Contents

Copyrights	1
Operator Safety	2
Warnings	3
Preface	5
Contents	6
Section I. Indications for Use	7
Section II. Contraindications and Warnings	7
Section III. Precautions	7
Section IV. Infection Control	9
Section V. Installation	12
Section VI. General Descriptions and Information of Parts.....	18
Section VII. Techniques	22
Section VIII. System Maintenance	24
Section IX. Troubleshooting	25
Section X. Specifications	26
Section XI. Disposal	26

SECTION I : INDICATIONS FOR USE

Ultrasonic procedures:

- Removal of calculus and plaque during dental prophylaxis.
- General supra and sub-gingival scaling applications.
- Periodontal debridement for all types of periodontal diseases.
- Endodontic procedures.

Polishing procedures:

- Removal of calculus and plaque during dental polishing.
- General supra and sub-gingival polishing applications.
- Periodontal debridement for all types of periodontal diseases.
- Endodontic procedures with optional bur tools.

Section II : Contraindications and Warnings

- Do not use the ART-SP6 for restorative dental procedures involving amalgam condensation.
- Do not use the ART-SP6 if the patient or operator has a pacemaker.
- Do not immerse the ART-SP6 in liquid. If the unit has fallen into liquid, do not reach into the liquid.
- Do not modify the ART-SP6. Modifications may violate safety codes and endanger the operator or subject, as well as void the warranty terms and conditions.

Section III : Precautions

Precautions for the ART-SP6 and all Ultrasonic Scaling systems


- Operation of the ART-SP6 must be performed by qualified veterinarians, technicians, and other trained personnel.
- Ensure that an adequate amount of water is provided during operation to cool and prevent the hand piece and piezo tip from overheating.
- Turn off the water source to the ART-SP6 when the scaler is not in use.
- Use caution to prevent transmission of contagious diseases as all ultrasonic scalers produce aerosols.
- Keep away from heat sources. Excessive heat may damage the electronic components.
- Seek advice from a physician prior to operating the ART-SP6 if the operator is pregnant.

Precautions for Ultrasonic Prophylaxis Procedures


- Discard and replace piezo tips as needed. Scaling tips wear out with use and should be replaced after 2mm of wear (approximately 3 to 4 months) to maintain optimal performance
- Discard and replace piezo tips that are bent, damaged, or reshaped. Do not reuse deformed tips
- BONART piezo tips are made of stainless steel. Avoid direct contact with teeth to prevent damage to the enamel. We strongly advise practicing on models and aluminum plates prior to operating on a live subject with the ART-SP6.
- Avoid direct contact of piezo tips with the subject's lips, cheeks, and tongue.
- Wipe accretions from the tooth using the long axis of the scaler tip. DO NOT gouge the tooth.

Precautions for hand-piece


Do Not Twist and Pull Cable


Do Not Bend


Do Not Twist

Do Not Bundle

Pull To Disconnect

SECTION IV : INFECTION CONTROL

Infection Control:

Always use standard personal protective equipment such as facemasks, eyewear, face shields, gloves, and protective gowns.

In addition to using protective equipment, ensure that your water supply system conforms to applicable Centers for Disease Control and Prevention (CDC) and American Dental Association (ADA) standards, and that all required flushing, chemical flushing, and general infection control procedures are followed.

To further protect against infection, using distilled or purified water is recommended.

Cleaning and Sterilization:

All sterilization procedures must be followed in accordance to the EN ISO17664 standards.

In order to avoid bacterial or viral infections, always clean and sterilize the following components after each treatment:

1. Hand piece
2. Tips
3. Metal Wrenches
4. Micro-motor
5. Straight hand piece (Straight nose cone)
6. Prophy head

The hand piece and tips, wrench, straight handpiece and prophy head can be autoclaved and can withstand a maximum temperature of 135°C or (275°F) for a **maximum of 3 minutes**. After completing the sterilization cycle, air-dry for about 20-30 minutes before reusing.

WARNING: *Do not use any other method of sterilization (dry heat, radiation, ethylene oxide, gas, low-temperature plasma, etc.)*

NOTE : The Sterilization Assurance Level (SAL) of steam autoclave should be 10⁻⁶ (According to ISO 13683: Sterilization of health Care products)

NOTE: Prior to cleaning and sterilizing your ART-SP6 unit, always remember to turn off the device by using the power switch and to disconnect the power plug from the outlet.

Cleaning the ART-SP6 unit

- Remove the hand piece from the cable, and the tip from the hand piece.
- Clean and disinfect the surface of the outer casing, cords and connectors with a cloth and a gentle disinfectant.

Cleaning and sterilizing the hand piece

- Remove the hand piece from the hand piece cable and the tip from the hand piece.
- Take special care not to break the threading pin of the hand piece. The threading pin is where the tip attaches to the hand piece.
- Disinfect the hand piece with a cloth and mild disinfectant containing a neutral pH.
- Dry the electric contacts by blowing air onto them with the syringe.
- Seal the hand piece in a sterilization pouch (without any tips). If you are short on sterilization pouches, you may use the surgical tray in place of the sterilization pouch.
- Autoclave the hand piece.

WARNING : Do not immerse the hand piece cable in sterilizing fluid. The ART-SP6 can be severely damaged if sterilization fluid gets inside the unit.

Cleaning Piezo Tip or Wrench

- Clean the tip or wrench preferably in an ultrasonic tank and rinse in distilled water.
- Dry the tip or wrench.
- Disinfect the tip with a cloth and mild disinfectant containing a neutral pH.
- Seal the tip or wrench inside a sterilization pouch. If you are autoclaving more than 1 tip, place each tip or wrench in an individual bag. You may also use the tip holder that comes with the ART-SP6 in place of the sterilization pouch.
- Autoclave the tip or wrench .

WARNING : Before starting the sterilization cycle, make sure that the inside of the tip or wrench is completely dry by blowing air through the internal hole with the syringe. This will prevent stains and patches from appearing on the surface of the tip or wrench.

Cleaning the micro-motor (optional accessory)

Remove all attachments to the micro-motor prior to cleaning. The outer surface of the micro-motor should be cleaned with an antiseptic solution. **Never oil the micro-motor. Ensure that no moisture gets inside the micro-motor.**

WARNING: Do not dip the micro-motor and cord directly into sterilization fluids. The chemicals will interfere with normal operations if it reaches the inside of the system.

WARNING: Do not leave disinfectants on parts longer than recommended by the manufacturer of the disinfectant to avoid damaging materials.

NOTE: The micro-motor should be cleaned after each use between patients.

Cleaning the straight hand piece (straight nose cone) attachment (optional accessory)

- Clean the straight hand piece preferably in an ultrasonic tank and rinse in distilled water.
- Dry the straight hand piece.
- Disinfect the straight hand piece with a cloth and mild disinfectant containing a neutral pH.
- Seal the straight hand piece inside a sterilization pouch. If you are autoclaving more than 1 tip, place each straight hand piece in an individual bag. You may also use the tip holder that comes with the KRUUSE-SP2 in place of the sterilization pouch.
- Autoclave the straight hand piece .

Cleaning the Propy Head (unless disposable) (optional accessory)

The prophy head should be cleaned to remove hair and debris approximately after every 10 dental polishing.

1. Remove prophy angle from handpiece and rubber cup.
2. Remove head gears by turning knurled knob clockwise (threads are reversed) turn to your right.
3. Clean with fine brush and lubricate head gears with mineral oil, Vaseline or any non-toxic lubricant.
4. Place head gears back and tighten by turning Left or counterclockwise.
5. Normal polishing speeds are around 2,000 to 3,000 RPM, exceeding this may damage your prophy angle.

SECTION V : INSTALLATION

General Information

If the installation of your ART-SP6 system is performed by non-Bonart dealer personnel, check to see that the requirements below are followed.

Water Line Requirements


- Incoming water supply line pressure to the Scaler must be 25 psi (172 kPa minimum) to 60 psi (414Kpa maximum). If your dental water supply line pressure is above 60 psi, install a water regulator on the water supply line to your ART-SP6 unit.
- After installing the water supply system, flush the water thoroughly prior to connecting to the scaler.
- Manually shut off the water valve of your dental water supply system prior to leaving the office unoccupied and unattended.

Electrical Requirements

Refer to Section X: Specifications

Unpacking the Unit

Carefully unpack your ART-SP6 unit and verify that all components and accessories are included:

(Figure 1) All components and accessories

Item	Quantity	Item	Quantity
(1) Main Unit	1	(2) Foot Control / Foot Switch	1
(3) AC Power Cord Set.	1	(4) PU Water Tubing	1
(5) User Manual & Literature Packet	1	(6) Straight Hand piece.	1
(7) Metal Torque-Limited Wrench	1	(8) ART-BS1 Universal Pointed Tip ART-BS2 Chisel Head Flat Tip ART-BS3 Round Flat Tip	2 1 1
(9) Micro-motor			Optional
(10) Micro-motor seat			Optional
(11) Straight Hand Piece			Optional

NOTE: Other piezo tips are available and can be purchased separately through dealers. Included items are subject to change depending on location.

NOTE : While unpacking, check your ART-SP6 unit for any damages. If any damages are found, please contact your dealer immediately. Enter the unit serial number on your warranty card and mail it within 14 days from the date of purchase.

Power Cord / Power Connection

- Always make sure that the power switch is set to the OFF position before performing the following tasks:
 - Plugging and unplugging the AC cord into the back of the unit.
 - Plugging and unplugging the 3-prong plug into a grounded electrical outlet.

Safety Instructions

A. Grounding:

Verify that the main unit is connected to a grounded outlet prior to connecting other parts to the main unit.

B. Main voltage range and fuse:

Check the voltage of the ART-SP6 for compatibility with your local voltage before plugging the power cord into the outlet.

Foot Control Cable Assembly Connection

Align the pins of the foot control plug with the receptacle on the back of the device and push firmly until seated.

Water Supply Line Connection

Push the blue water tube (hose) over the stainless steel receptacle on the back of the unit. Then tighten the screw over it.

Connect the quick connector to the water supply line. The ART-SP6 standard package comes with a male quick connector. If a female quick connector is necessary, please contact a BONART authorized dealer to purchase it.

Inspect all connections for leaks


To remove the water line from the ART-SP6, turn off the water supply or disconnect the water supply line from the water tube. Then loosen the screw on the water tube from the receptacle of the unit and gently pull the water tube off.

Hand piece Cable Assembly Connection

The hand piece is the housing for the piezo tips. It is the part that you hold during operation. The hand piece must be attached to the hand piece cable prior to using. To attach, connect the free end of the cable to the matching receptacle on the back end of the hand piece. Do not twist or turn.

To detach, pull the cable from the hand piece in a straight motion. Do not twist, turn, or yank the hand piece from the cable.


Polishing Hand Piece Assembly and Start Up.

(Figure 2)

To assemble the polishing hand piece, connect the Prophy Angle, Straight Hand Piece, and Micro-motor by following the steps below:


1. Connect the straight hand piece and the micro-motor by pushing the straight hand piece over the rod of the micro-motor. It should snap into place.

(Figure 3)


2. Attach the prophy angle to the straight hand piece by following the steps below:

- A) Unlock the prophy angle from the straight hand piece by turning the locking ring clockwise. This should open the lock. Note: When you first unpack the straight hand piece, there may be a metal rod at the top that serves as a placeholder. Remove the rod after you open the lock.


(Figure 4)

- B) Match the groove at the base of the prophy head with the top of the straight hand piece and connect them together.

(Figure 5)

- C) Close the lock by turning the locking ring counter-clockwise. The prophy head should now be locked onto the straight hand piece and should not fall off during operation.

(Figure 6)

- D) If you wish to remove the prophy head, follow step A to release the lock on the prophy head.

WARNING : The prophy head should be cleaned and lubricated on a regular basis.

To polish, follow the steps below:

1. Place a rubber cup on the prophy head and put paste in the cup to begin polishing.
2. Press the Micro-motor (mode) button to switch from scaling to polishing mode
3. Adjust the Oscillating dial to a low speed, 2,000 rpm.
4. Press on the foot switch to engage the polishing hand piece and begin polishing.
5. Add paste as needed.

Assembling the Scaler Unit

- a. Ensure that the power switch is set to off. If plugged in, the power indicator light should not be lit.
- b. Plug the power cord of the main unit into a grounded AC power outlet.
- c. Insert the tip into the hand piece and lock it into place by using the enclosed plastic wrench. Be careful not to overturn the tip as it may cause the tip and hand piece to break.
- d. Ensure that the water tube connection to the water supply source and main unit is secure.

SECTION VI : GENERAL DESCRIPTION AND INFORMATION OF PARTS


Figure 7. ART-SP6 diagram

1. Main Unit

The main unit generates power and produces a signal that is passed to the hand piece. The hand piece is then powered up and vibrates the installed tip.

2. Hand piece

The hand piece is the housing for the piezo tip. It consists of four piezoelectric plates that vibrate. The plates are driven by the electrical signal from the Main unit. The ART-SP6 unit is compatible with Bonart BS-type piezo tips and other brands of piezo tips with a comparable thread. Do not use Bonart BE-type piezo tips with the ART-SP6.

3. Mode/Power control knob

The Mode/Power Control Knob is used to switch between polishing and scaling mode and control the power in scaling and/or polishing mode.

Do not press the foot switch while turning on the ART-SP6

4. Micro-motor Socket

The micro-motor socket is where you connect the micro-motor to the main unit. It is located in the front of the main unit at the bottom left.

5. Digital Display

The digital display shows the current operation mode (polishing or scaling mode) and the power / micro-motor speed level.

6. Water control knob

The water control knob allows you to control the amount of water that flows through the hand piece. If turned clockwise, the amount of water will decrease. If turned counter-clockwise, the amount of water will increase. Please note that water is necessary to cool the tip and prevent overheating so it is necessary to ensure that an adequate amount of water flows through the hand piece while operating the device. The larger the water stream, the lower the temperature, and vice versa.

WARNING : Do not turn the water knob counter-clockwise more than three full turns. Doing so will damage the internal water line needle.

WARNING : If needed, gently adjust the water control knob to prevent any serious damages to the water valve.

7. Water hose input connector

The water hose input connector is the contact point for the water tube and the main water supply line. Use either a male or female quick connect to connect the water tube to the water supply.

8. Footswitch Socket

The foot switch socket is where you connect the foot switch to the main unit. It is located next to the water hose input on the back of the main unit.

9. Power cord socket

The power cord socket is where you connect the power cord to the main unit.

10. Power button (ON/OFF)

The power button is used to turn the ART-SP6 main unit on and off. When the power is on, the LCD indicator should be lit. When turned off, the LCD power indicator should not be lit.

Do not press the foot switch while turning on the ART-SP6.

11. Forward / Reverse switch

The Forward/Reverse switch allows you to control the direction of the micro-motor. It is located on the back of the main unit.

The function of each button or knob on the control panel is detailed below.


Figure 8. ART-SP6 Front-panel


1. Mode/Power control knob

The Mode/Power Control Knob is used to switch between polishing and scaling mode and control the power in scaling and/or polishing mode. When the power is switched on, the unit will start in scaling mode. Push the knob to switch polisher mode. In scaling mode, turn the knob to adjust the scaling power from 1 to 50. In polishing mode, turn the knob to adjust the micro-motor speed from 30 to 300 (x100 R.P. M.) .


2. Digital Display

When the power is switched on, you should see a 10 second countdown on the screen as the system loads scaling mode. The most recent settings will also be displayed on the screen.

While you are in scaling mode, the border of the Scaler icon will be bolder, and the numbers below will display your current power level.

When you are in polishing mode, the border of the Polisher icon will be bolder, and the numbers below will display your current micro-motor speed (x100 R.P.M.).

3. Water control knob

The water control knob allows you to control the amount of water that flows through the hand piece. If turned clockwise, the amount of water will decrease. If turned counter-clockwise, the amount of water will increase. Please note that water is necessary to cool the tip and prevent overheating so it is necessary to ensure that an adequate amount of water flows through the hand piece while operating the device. The larger the water stream, the lower the temperature, and vice versa.

4. Micro-motor Socket Seat

The micro-motor connects to the main unit through this socket.

SECTION VII: TECHNIQUES

Performing Ultrasonic Scaling Procedures

1. Use purified or distilled water to prevent infection when patients experience tissue laceration during treatment
2. Keep the power cord tidy to avoid tripping and other accidents.
3. Position the footswitch in an easily accessible spot for the user. Keep the footswitch cord tidy to avoid tripping and other accidents.
4. Lock the tip in the hand piece securely using the enclosed wrench.

CAUTION: DO NOT FORCE OR OVERTURN THE TIP INTO PLACE! Doing so will damage the tip, hand piece, and wrench.

NOTE: the scaler will not operate if the tip is either removed from the hand piece or improperly installed.

5. Hold the hand piece over a sink or drain with a tip installed. Verify that water is reaching the tip.
6. Check your Bonart piezo tip for wear and replace as needed.
7. Use the foot control to regulate water flow to the hand piece. Hold down the foot control to allow water to flow through. Release the foot control to stop.
8. Adjust the water flow to the desired level using the water-adjustment knob. Increase the water flow by turning the water control knob clockwise. Decrease the water flow by turning the water control knob counter-clockwise.





NOTE: Increased water flow results in a cooler hand piece temperature and vice versa.

CAUTION: During operation, a continuous flow of water is required to keep the hand piece cool.

Applications of the TIP





A. Tangential application (BS1 TIP)

Do not apply the BS1 tip directly to the tooth. Doing so may damage the enamel. Control the hand piece with a slow and steady motion.





B. Frontal application (BS2 TIP)

Apply the BS2 tip directly against the tartar, but not directly against the tooth. Use a slight amount of pressure.

C. Tangential application (BS3 TIP)

Do not apply pressure when applying the BS3 tip to the enamel.

SECTION VIII: SYSTEM MAINTENANCE

Hand piece maintenance

1. After each use the hand piece and tip should be rinsed with clean water for about 20 to 30 seconds.
2. Inspect the hand piece cable daily to ensure it is in good condition.
3. For sterilization of the hand piece please refer to Section IV.

TIP maintenance

Worn out tips can adversely affect performance, resulting in insufficient power and vibration. Check tips regularly for wear and tear and replace as necessary.

Main Unit Maintenance

The housing of the ART-SP6, as well as the hand piece cord, should be cleaned and disinfected (using alcohol or soap and warm water) on a daily basis.

SECTION IX: TROUBLESHOOTING

While Bonart offers repair service, the following trouble shooting procedures can help you save time.

Troubleshooting Guide

PROBLEM	POSSIBLE CAUSE	SOLUTION
The device does not turn on when the switch is in the "ON" position.	The power cord is not plugged in. The unit may be defective.	Check your connection to the outlet and the unit.
The device is turned on, but does not vibrate when you press the foot switch.	Tip is worn-out or is not secured properly.	Check that the tip is not worn-out and ensure that it has been properly secured to the hand piece with the included tip wrench. Replace the tip if old or damaged.
	The foot switch may be broken or defective.	Reposition the foot switch and cable. If functionality is intermittent, a short in the wire may be present. Replace the foot switch if necessary.
The device is turned on, but does not produce water when you press the foot switch.	The water supply is not connected to the unit properly.	Check your water supply for sufficient flow, and ensure that the water tube is properly secured to the back of the unit.
	The filter may be clogged.	Clean or replace the filter
Insufficient vibration/oscillation.	The tip may be old, worn-out, or not properly inserted.	Check the tip for wear and tear, and ensure that it is properly inserted by using the included tip wrench. Replace the tip if it is old or damaged.
	The power setting is too low.	Turn the power dial to a higher setting.
Insufficient water pressure.	The filter may be clogged.	Replace the filter screen or the filter itself.
	The pressure from the water supply source may be too low.	Check that pressure of the water supply source is between 25 and 60 PSI.

SECTION X: SPECIFICATIONS

The ART-SP6 scaler is designed and manufactured to meet the most demanding environment. It's specifications are listed below

Classification

- Protection against electric shock : Class I
- Degree of protection against electric shock : Type BF
- Mode of operation : Continuous
- Medical device directive : IIa

Standards Followed

ISO13485: Medical Devices, Quality Management Systems, Requirements for regulatory purposes

93/42/EEC: MDD CE Directive

EN60601-1: Medical Electrical Equipment, General Requirement for Safety

EN60601-1-2: Medical Electrical Equipment, Electromagnetic Compatibility, Requirements and Tests

● Specification

- Power supply
115V \pm 5% ~50/60Hz 125VA
230V \pm 5% ~50/60Hz 125VA
- Working frequency
26KHz ~ 32KHz
- Polisher speed
3000~30000 R.P.M
- Water supply
25~60 PSIG(172~414KPa)
- Dimension
23cm(L) x 21cm(W) x 9cm(H)
- Weight
3.7 Kg (including the handpiece)
- Hand piece Cable
250 cm
- Footswitch Cable
250 cm

Operation environment

- Temperature
0 ~ 40 °C
- Relative Humidity
10% ~ 90% (non-condensing)

Transport and storage conditions

- Temperature
0 °C ~ 60°C
- Relative Humidity
10% ~ 90% (non-condensing)
- Atmospheric pressure
860~1060 hPa

SECTION XI: DISPOSAL

Please follow county and state regulations for disposal of the ART-SP6.

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COPYRIGHT© 2014 BONART CO. , LTD.