iM3 Bioceramic RCS

Canal Filling and Sealing Material

For veterinary use only

A ready to use injectable premixed filling and sealing material which is based on bioceramic technology. It is an alumina free calcium silicate based material with excellent handling characteristics and superior physical properties. It is a radiopaque material which requires water to set and harden.

Composition:

Calcium Silicates, Calcium Phosphates, Zirconium Oxide, Bioactive Glasses, Calcium Sulphate, Calcium Oxide, Fillers, Accelerators and thickening agents.

Indications for use:

1. Permanent obturation /filling of root canals.
2. As a root canal sealer for teeth of the secondary dentition in combination with root canal points. (suitable for single cone and lateral condensation).

Directions for use:

The cement requires no mixing as it is supplied as a premixed paste in a preloaded syringe with disposable intracanal delivery tips. Prior to the application of the material, the prepared root canal should be dried using a paper point. Avoid over drying as the sealer needs moisture to completely set and harden.

1. Sealer technique:

A. Conventional method:

The sealer can be dispensed from the syringe on a mixing pad and delivered to the root canal space by carrying it with a gutta-purca point or an endodontic instrument.

B. Tip delivery technique:

Remove the syringe cap and attach the intracanal delivery tip. Insert the tip into the coronal third of the root canal space and smoothly dispense a small amount of material by compressing the syringe plunger. Coat the walls of the canals evenly using a #15 endodontic file or a gutta-purca point. In case of lateral condensation place additional gutta-purca points coated with the sealer. After obturation the gutta-purca is cut at the orifice using a heat source. Excess sealer in the root canal chamber should be wiped off using a cotton balls or swabs.

2. Filling technique:

The material can also be used for filling the root canal space independently. In that case, proper care should be taken to prevent extrusion of material. An apical limiting barrier and is strongly advised in such cases.

Contraindications:

Hypersensitivity against calcium silicates or other components of the root canal filling material.

Warnings:

A. The sealer which may cause sensitization in susceptible patients. Do not use in patients allergic to calcium silicates or any other components of this product.

B. Avoid eye contact to prevent irritation. In case of contact with eyes, rinse immediately with plenty of water and seek

medical attention. Avoid skin contact to prevent possible allergic response. Unset material may cause irritation and therefore contact with skin and mucous membrane should be avoided. In case of contact, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.

C. Do not swallow or take internally. If accidental swallowing occurs, drink lot of water. If nausea or illness develops, seek medical attention.

Avoid overfilling. When overfilling occurs, the material is usually tolerated very well by the surrounding tissue. The material should be immediately removed in cases of extrusion in to the mandibular canal.

Storage conditions

This product should be used at room temperature. Contact with water should be strictly avoided. Inadequate storage conditions could shorten shelf life and product performance. We recommend to store at temperature between 10°C (50°F) and 24 °C (75°F) and keep syringe cap tightly closed.

Shelf Life - 3 years from the date of manufacturing.

Disposal conditions:

To be disposed off based on the prevailing rules of the country in which it is used.

USA (800) 664 6348 info@im3usa.com
Australia +612 9420 5766 sales@im3vet.com
Ireland +353 (0)1 6911277 sales@im3dental.com

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